

AD-A116 804

ARIZONA UNIV TUCSON DEPT OF CHEMISTRY  
PAIRED ION CHROMATOGRAPHY OF CERTAIN MONOSACCHARIDES.(U)  
SEP 81 T GNANASAMBANDAN, H FREISER  
TR-5

F/6 7/3

N00014-81-K-0576

UNCLASSIFIED

NL

[RE]  
AL-1  
01104



END  
DATE  
FILMED  
6-82  
DTIC

REPORT DOCUMENTATION PAGE		1. READ INSTRUCTIONS BEFORE COMPLETING FORM
1. REPORT NUMBER 5	2. GOVT ACCESSION NO. AD-A116804	3. RECIPIENT'S CATALOG NUMBER
4. TITLE (and Subtitle) "Paired Ion Chromatography of Certain Mono-saccharides"		5. TYPE OF REPORT & PERIOD COVERED Technical -- Interim
7. AUTHOR(s) T. Gnanasambandan and H. Freiser		6. PERFORMING ORG. REPORT NUMBER
9. PERFORMING ORGANIZATION NAME AND ADDRESS University of Arizona Department of Chemistry Tucson, Arizona 85721		8. CONTRACT OR GRANT NUMBER(s) N00014-81-K-0576
11. CONTROLLING OFFICE NAME AND ADDRESS Materials Sciences Division Office of Naval Research Arlington, Virginia 22217		10. PROGRAM ELEMENT, PROJECT, TASK AREA & WORK UNIT NUMBERS
14. MONITORING AGENCY NAME & ADDRESS (if different from Controlling Office)		12. REPORT DATE
		13. NUMBER OF PAGES 5
		15. SECURITY CLASS. (of this report) unclassified
		15a. DECLASSIFICATION/DOWNGRADING SCHEDULE
16. DISTRIBUTION STATEMENT (of this Report) This document has been approved for public release and sale; its distribution is unlimited		
17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)		
18. SUPPLEMENTARY NOTES to be published in Analytical Chemistry		
19. KEY WORDS (Continue on reverse side if necessary and identify by block number) monosaccharide chromatographic analyses, fructose, ribose, and xylose		
20. ABSTRACT (Continue on reverse side if necessary and identify by block number) A paired ion dye-impregnated reversed phase chromatographic column is shown to give well defined chromatographic zones for monosaccharides as their borate complexes, permitting determination of sub-microgram levels of fructose, xylose, and ribose.		

DD FORM 1473

JAN 73

EDITION OF 1 NOV 65 IS OBSOLETE  
S/N 0102-LF-014-6601

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

AD A116804

DTIC FILE COPY

DTIC  
ELECTE  
JUL 12 1982

D

H

OFFICE OF NAVAL RESEARCH  
Contract N00014-81-K-0576

TECHNICAL REPORT NO. 5

"Paired Ion Chromatography  
of Certain Monosaccharides"

by

T. Gnanasambandan and H. Freiser

submitted to

Analytical Chemistry

University of Arizona  
Department of Chemistry  
Tucson, Arizona 85721

Reproduction in whole or in part is permitted for  
any purpose of the United States Government

This document has been approved for public release  
and sale; its distribution is unlimited

## Paired ion Chromatography of Certain Monosaccharides

Sir:

Several chromatographic methods (1,2,3) have been reported for separation of sugars but this still remains a formidable task for determination of trace levels. The packings which have been most successful in separating sugars are those with amine groups bound to inert silica or synthetic supports such as U-Bondapak carbohydrate (Waters Assoc.), Lichrosorb-NH<sub>2</sub> (Merck) and Aminex resins (Bio-Rad Lab) (4). These methods involve elution with water/acetonitrile system or water at 80-85°C.

In the course of our study in reverse phase ion-pair HPLC (5) we found that sugars such as fructose, can be ionized on column and paired with a chromophore to effect improved separation and detection.

One of the interesting aspects of this separation is that the non-ionic species react on the column to form suitable species and separated. Generally the analyte is subject to chemical pre-treatment before separation. This is the first report of an on-column reaction in reverse phase HPLC.



Accession For	
NTIS	CRUI
DTIC	TR
Unannounced	
Justification	
By	
Distribution	
Availability Codes	
Dist	Special

A

## Experimental

A modular chromatographic unit consisting of an Altex Pump (Model 100A), a variable-wave length uv-vis detector (Schoeffel Instruments No. 770) and a Fisher Recordall Series 500 strip chart recorder was used. Absorbance measurements were made at 651nm. A 25cm Zorbax ODS column, described by the manufacturer as having 15% carbon loading, was employed.

The mobile phase used was methanol-(0.1M aqueous borate buffer PH 7.5) in 2:98 V/V% containing  $1 \times 10^{-4}$  M methylene blue (as chloride).

Before the chromatographic separations were undertaken, the column was conditioned by passing sufficient mobile phase to "saturate" it with methylene blue (6). The flow rate used for both column saturation and subsequent chromatography was 1.0ml per min. (1% W/V) solutions of the fructose, xylose and ribose (separately) in the mobile phase were added through 10 $\mu$ l loop sampling valve (rotary valve injector SP-419-0410). As may be seen from the data, reasonable separation of the three sugars tested was obtained at sub-microgram detection limits ( $\sim 10^{-7}$  g).

Table 1. Capacity Factors of Monosaccharides

Sugars	k'
Fructose	3.5
Ribose	2.9
Xylose	1.2
Mobile phase:	0.1M Borate Buffer pH7.5 Methylene blue $10^{-4}$ M 2% MeOH/H <sub>2</sub> O V/V
Flow rate:	1.0ml per min.
Sample:	10 $\mu$ l 1.0% wt/V Solutions

The use of borate complex formation to produce anionic sugar species results in a dramatic improvement in selectivity. The formation constant for the glucose-borate complex is far smaller ( $\sim 10^2$ ) than the corresponding fructose complex ( $\sim 10^4$ ) (7). As a result, glucose does not give any detectable response under the experimental conditions employed here, making the determination of fructose entirely free from glucose interference. Changing these conditions to accomodate the lower stability of the glucose complex, e.g. by increasing borate ion concentration, would very likely permit determination of glucose also.

Although further work is required to clearly establish the optimum conditons necessary to separate all sugars, it would appear that reactions occur on the column to give suitable products for effective separation and detection. Further work along these lines is under way in this laboratory.

#### Literature cited

1. Ladish, M.R., Tsao, G.T., J. Chromatogr., 1978, 166, 18-100.
2. Palla, G., Anal. Chem., 1981, 53, 1966-1967.
3. Linden, J.L., Lawhead, C.L., J. Chromatogr. 1975, 105, 125-233.
4. Conrad, E.C., Palmer, J.K., Food Technol (Chicago), 1976, 84-92.

5. Gnanasambandan T., Freiser, H., Anal. Chem., 1981, 53, 909-911.
6. Gnanasambandan, T., Freiser, H., Anal. Chem., in press.
7. Roy, G.L., Laferriere and Edwards, J.O., J. Inorganic Nuclear Chem., 1957, 4, 106.

Acknowledgment. This work was supported by the Office of Naval Research.

T. Gnanasambandan

H. Freiser

Department of Chemistry  
University of Arizona  
Tucson, Arizona  
April 19, 1982



MAILING LIST FOR OFFICE OF NAVAL RESEARCH REPORTS

Mr. John Boyle  
Materials Branch  
Nabal Ship Engineering Center  
Philadelphia, PA 19112

Dr. H. Chernciff  
Department of Mathematics  
Massachusetts Institute of Technology  
Cambridge, MA 02139

Walter G. Cox, Cod 3632  
Naval Underwater Systems Center  
Building 148  
New Port, RI 02840

Defense Technical Information Center (12)  
Building 5, Cameron Station  
Alexandria, VA 22314

Dr. M.B. Denton  
Department of Chemistry  
University of Arizona  
Tucson, AZ 85721

The Assistant Secretary of the Navy  
(RExS)  
Department of The Navy  
Room 4E736, Pentagon  
Washington, DC 20350

Department of Physics & Chemistry  
Naval Postgraduate School  
Monterey, CA 93940

Dr. John Duffin  
United States Naval Postgraduate School  
Monterey, CA 93940

Dr. Kent Eisentraut, MBT  
Air Force Materials Laboratory  
Wright-Patterson AFB, OH 45433

Dr. Christie G. Enke  
Michigan State University  
Department of Chemistry  
East Lansing, MI 48824

Dr. G.M. Hieftje  
Department of Chemistry  
Indiana University  
Bloomington, IN 47401

Mr. James Kelley  
DTNSRDC Code 2803  
Annapolis, MD 21402

Dr. Bruce Kowalski  
Chemistry Department  
University of Washington-Seattle  
Seattle, WA 98105

Dr. Rudolph J. Marcus  
Office of Nabal Research  
Scientific Liaison Group  
American Embassy  
APO San Francisco, CA 96503

Dr. George H. Morrison  
Department of Chemistry  
Cornell University  
Ithaca, NY 14853

Naval Civil Engineering Laboratory  
Attn: Dr. R.W. Drisko  
Port Hueneme, CA 93401

Naval Weapons Center  
Attn: Dr. A.B. Amster  
Chemistry Division  
China Lake, CA 93555

Naval Ocean Systems Center  
Attn: Dr. S. Yamamoto  
Marine Sciences Division  
San Diego, CA 91232

Naval Ocean Systems Center  
Attn: Mr. Joe McCartney  
San Diego, CA 92152

Commander, Naval Air Systems Command  
Attn: Code 310C (H. Rosenwasser)  
Department of the Navy  
Washington, DC 20390

Director, Naval Research Laboratory  
Attn: Code 6100  
Washington, DC 20390

Naval Ship Research & Development Center  
Attn: Dr. G. Bosmajian  
Applied Chemistry Division  
Annapolis, MD 21401

Office of Naval Research (2)  
Attn: Dr. Richard S. Miller  
800 N. Quincy Street  
Arlington, VA 22217

Office of Naval Research  
Attn: Code 472  
800 N. Quincy Street  
Arlington, VA 22217

ONR Branch Office  
Attn: Dr. George Sandoz  
536 S. Clark Street  
Chicago, IL 60605

ONR Eastern/Central Regional Office  
Attn: Dr. L. H. Peebles  
Building 114, Section D  
666 Summer Street  
Boston, MA 02210

ONR Western Regional Office  
1030 East Green Street  
Pasadena, CA 91106

Dr. R. A. Osteryoung  
Department of Chemistry  
SUNY, Buffalo  
Buffalo, NY 14214

Dr. S. P. Perone  
Department of Chemistry  
Purdue University  
West Lafayette, IN 47907

Dr. Victor L. Rehn  
Naval Weapons Center  
Code 3813  
China Lake, CA 93555

Dr. Freed Saalfeld  
Chemistry Division, Code 6100  
Naval Research Laboratory  
Washington, DC 20375

U.S. Army Research Office  
Attn: CRD-AA-IP  
P. O. Box 1211  
Research Triangle Park, NC 27709

Dr. A. L. Slafkosky  
Scientific Advisor  
Commandant of the Marine Corps  
(Code RD-1)  
Washington, DC 20380

Dr. D. L. Venezky  
Naval Research Laboratory  
Code 6130  
Washington, DC 20375

Dr. Isiah H. Warner  
Texas A&M University  
Department of Chemistry  
College Station, TX 77840

Dr. K. Wilson  
Department of Chemistry  
University of California, San Diego  
La Jolla, CA 92037

Dr. A. Zirino  
Naval Underseas Center  
San Diego, CA 92132

James R. Delaney  
Contract Specialist  
Office of Naval Research  
422 University District Bldg.  
1107 NE 45th Street JD-27  
Seattle, WA

9/23/81